TECHNICAL DATA SHEET



Dental tray filter paper



MATERIAL	 Cellulose 55 g/m² Chlorine and optical brightening agents free
PRODUCT INFORMATION	 Gas and steam absorbent, it avoids any contact between the instruments and the dental tray for a perfect instruments sterilization and preservation Colour: white Tear index, dry machine direction: >25 N Wet strength machine direction: >5.5 N Tensile stretch, dry: >6.0% General tolerance: +/- 10 %
PACKAGING	 Packed in polyethylene bag to protect from dust and humidity Pack sizes: 275 x 175 x 45 mm, 0.75 kg Export boxe sizes: 365 x 290 x 235 mm Should be store in a clean and dry place, away from sunlight

STANDARDS PACKAGING

Reference	Dimensions	Pieces/Pack	Packs/Box	Boxes/Pallet
TRAYW	18 x 28 cm	250	10	80