## TECHNICAL DATA SHEET



## Napkins



MATERIAL	<ul> <li>1 layer cellulose wadding (25-28gr/m²)</li> <li>1 layer polyethylene (10-12mμ)</li> <li>Polyethylene thermo paired, embossed and folded into napkins for better absorption of liquids</li> <li>Solvent free</li> </ul>
PRODUCT INFORMATION	<ul> <li>Standards utilization: to protect the patient</li> <li>Waterproof polyethylene</li> <li>Absorbent paper (hygienic)</li> <li>Soft and comfortable</li> <li>Shelf life: 10 years from the date of manufacture</li> <li>Class I medical device, single use product</li> <li>General tolerance: +/- 10 %</li> </ul>
PACKAGING	<ul> <li>Packed and protected from dust and humidity</li> <li>Different colours on request and depending the quantity</li> <li>Export box dimensions: 35 x 25 x 25 cm</li> <li>Should be store in a clean and dry place, away from the sunlight</li> </ul>

## **STANDARDS PACKAGING**

Reference	Description	Colour	Pieces/Bag	Bags/Box	Boxes/Pallet
NAPW	Napkins 45x33 cm	White	125	4	80
NAPB	Napkins 45x33 cm	Blue	125	4	80