## TECHNICAL DATA SHEET




## Alan-Silk suture




MATERIAL	<ul> <li>Composed of an organic protein (Fibroin)</li> <li>Black braided silk</li> <li>Needle made of stainless steel, silicone coated</li> </ul>		
PRODUCT INFORMATION	<ul> <li>Sterile surgical suture</li> <li>Non absorbable</li> <li>Multifilament</li> <li>Needle eases the passage without damaging the tissue</li> <li>Shelf life: 5 years from production date</li> <li>Class III medical device sterime, single use product</li> <li>General tolerance: +/- 10 %</li> </ul>		
PACKAGING	<ul> <li>12 sterile pouch/box to protect from dust</li> <li>Alan brand only</li> <li>Box dimensions: 62x130x62 mm</li> <li>Export box dimensions: 60x40x60 cm, 216 boites/carton</li> <li>Should be stored in a clean and dry place, out of direct sunlight and below a temperature of 30°C</li> </ul>		

## **STANDARDS PACKAGING**

Reference	Needle profile	Needle description	Needle length	Size	Thread length
SUTSI3	3/8	Reverse cutting	19 mm	3/0	70 cm
SUTSI31/2	1/2	Reverse cutting	17.4 mm	3/0	45 cm
SUTSI4	3/8	Reverse cutting	19 mm	4/0	70 cm
SUTSI41/2	1/2	Reverse cutting	17.4 mm	4/0	45 cm







## **ALAN & CO**

Avenue André Ernst 17 – 4800 Verviers – Belgium Tel : +32 87 31 10 59 – Fax : +32 87 31 61 69